

# HL-dsDNase

Best Choice for  
dsDNA Removal from  
RNA & Protein Samples!

NEW  
PRODUCT

Ideal for procedures where dsDNA presence is undesirable:

- RNA and protein samples rapid purification,
- PCR, qPCR Master Mixes, or other diagnostic kits decontamination.

- Low-temperature activity to protect your RNA or proteins!
- Irreversible inactivation in low temperature (52°C) to secure biological products' integrity!

## Overview

**HL-dsDNase** is a 43.3 kDa heat-labile endonuclease, originates from a cold water eukaryotic organism, recombinantly expressed in *Pichia pastoris*. The enzyme displays high specific activity solely towards double-stranded DNA leaving single-stranded DNA or RNA undamaged.

**HL-dsDNase** characterizes high specific activity and it is easily inactivated by heat treatment in moderate temperatures. It is intended for applications where the presence of dsDNA influences experiments' results in thermo-sensitive applications and is extremely useful for a rapid and safe digestion of genomic DNA in samples containing RNA or recombinant proteins.

## Applications

**Digestion of dsDNA (plasmid DNA, genomic DNA, etc.),**

- RNA and protein samples rapid purification,
- PCR, qPCR Master Mixes, or other diagnostic kits decontamination.

## Features and advantages

- **Highly active in the broad temperature range (optimum at 10-47°C),**
- **Highly active in typical buffer formulations and broad pH range (optimum at 7-8),**
- **Inactivation at moderate temperature (15 min. at 52°C, 1 mM DTT),**
- **The activity towards dsDNA is minimum 1000 times higher than towards ssDNA or RNA.**

Variable / Parameter	Activity Range	Optimum Activity
pH	6-10	7-8
Temperature	4-47°C	10-47°C
Mg <sup>2+</sup>	0-50 mM (ions Ca <sup>2+</sup> increase the activity)	0-20 mM
Ammonium sulfate	0-100 mM	0-50 mM
NaCl/KCl	0-250 mM	0-100 mM
Imidazole	up to 400 mM	0-300 mM
Urea	up to 400 mM	0-2 M
Glycerol	up to 50%	0-40%
Triton X-100	up to 2%	0-2%
DTT	0-100 mM	0-100 mM
β-merkaptotanol	0-2.5%	0-1%