

Tth DNA Ligase



1.2020



blirt

Tth DNA Ligase

Tth DNA Ligase catalyzes the NAD-dependent formation of phosphodiester bonds between adjacent 3' hydroxyl and 5'-phosphate termini in double stranded DNA. It is not active against single stranded DNA or RNA and blunt ended DNA. Enzyme is isolated from *Escherichia coli* strain containing plasmid carrying the *Thermus thermophilus* DNA ligase gene. **Tth DNA Ligase** is stable and active in optimum ligation temperature range of 45–65°C, which is 7–10°C higher than that of T4 DNA ligase. The final reaction ligation temperature is determined by the T_m of the substrates. High ligation temperature eliminates the nonspecific ligation.

Features and advantages

- Stable and active in high temperatures
- Highly specific ligation
- LCR (Ligase Chain Reaction)
- LDR (Ligase Detection Reaction)
- NGS (Next-Generation DNA Sequencing)
- RED (Repeat Expansion Detection)
- RCA (Rolling Circle Amplification)
- PLA (Proximity Ligation Assay)

Note: Some applications in which this product may be used may be covered by patents or patent applications applicable in certain countries. Because purchase of this product does not include a license to perform any patented application, users of this product may be required to obtain a license depending upon the particular application and country in which the product is used.

Stability

Enzyme retains full activity after incubation for 1 week at 37°C. The half-life of enzyme is about 48 hours at 65°C. 10x Tth Ligation Buffer is stable for 1 week at 37°C. Up to twenty freeze/thaw cycles will not compromise 10x Tth Ligation Buffer performance.

Protocol

1. Add the reaction reagents listed below to a sterile nuclease-free tube. The reaction agents should be added in the following order:

Component	Volume
Nuclease-free water	up to 25 μ l
10x Tth Ligation Buffer	2.5 μ l
Tth DNA Ligase 5 U*/ μ l (75 CEU/ μ l)	0.5–1 μ l
DNA	0.5–1 μ g

2. Mix gently and spin briefly.
3. Incubate for 10 min at 45–65°C. Optimum ligation temperature range is determined by the T_m of the substrates.

Quality control

Tth DNA Ligase activity is tested in reaction with bacteriophage lambda DNA digested with Sall and SmaI, with a dilution series of ligase. Results are assayed by agarose gel electrophoresis. Free of unspecific DNA and RNA nucleases contamination.

Storage Buffer

50 mM Tris-HCl (pH 7.5), 100 mM NaCl, 0.1 mM EDTA, 0.1% Triton X-100, 1 mM DTT, 50% glycerol

10x Tth Ligation Buffer

200 mM Tris-HCl (pH 8.3), 250 mM KCl, 100 mM MgCl₂, 5 mM NAD, 0.1% Triton X-100

Tth DNA Ligase

Components	EN13-025 250 U (3750 CEU)	EN13-250 2500 U (37 500 CEU)	EN13-S 25 U (375 CEU)
Tth DNA Ligase 5 U*/ μ l (75 CEU/ μ l)	50 μ l	500 μ l	5 μ l
10x Tth Ligation Buffer	125 μ l	1250 μ l	12.5 μ l

* Unit definition may vary between manufacturers

Additional information

Unit definition

One unit of Tth DNA Ligase catalyzes the ligation of 50% of the cos sites present in 1 μ g of bacteriophage lambda DNA in 1 minute at 45°C.

1 U (Unit) of Tth DNA Ligase = **1 Ampligase® Unit** = **15** cohesive end units (CEU).

Storage conditions


All components should be stored at -20°C.

Up to twenty freeze/thaw cycles will not compromise 10x Tth Ligation Buffer performance.

Shipping conditions

Shipping on dry or blue ice.

Ampligase® is a registered trademark of Epicentre

 For research use only

Expiry

Information on the label.