

ATTACHMENT TO

IX. SAMPLE PREPARATION

Below you will find additional information to sample preparation which should ensure high isolation efficiency and pDNA purity.

EXTRACTME PLASMID MAXI Kit (EM18)

EXTRACTME PLASMID MAXI ENDOTOXIN-FREE Kit (EM19)

		Recommended Culture Volume according to a culture density A_{600}		
Plasmid Type	Pellet Wet Weight	OD600=2	OD600=4	OD600=6
High-copy number	3 g	600 ml	300 ml	200 ml
Low-copy number	6 g	1000 ml	600 ml	400 ml