


DNA & RNA Extraction Kits

*Fast and efficient isolation of nucleic acids
from a wide range of biological materials*



MICRO SPIN KITS




**TOTAL RNA
MICRO SPIN**
EM31




**DNA CLEAN-UP &
GEL-OUT MICRO SPIN**
EM28




MINI SPIN KITS




RNA & DNA
EM15



TOTAL RNA
EM09.1




TOTAL RNA PLUS
EM11.1




miRNA
EM12




RNA BACTERIA & YEAST
EM25




DNA TISSUE
EM03



DNA TISSUE PLUS
EM04



GENOMIC DNA
EM13



DNA BLOOD
EM05



DNA SWAB & SEMEN
EM06



PLASMID MINI
EM01.1



DNA BACTERIA
EM02




DNA CLEAN-UP
EM07.1




DNA GEL-OUT
EM08.1




**DNA CLEAN-UP
& GEL-OUT**
EM26.1



**MIDI / MAXI
COLUMNS**



**PLASMID MIDI /
ENDOTOXIN-FREE**
EM16/17



**PLASMID MAXI /
ENDOTOXIN-FREE**
EM18/19

USE OUR TABLE to find a perfect match for your extraction procedures or ASK OUR EXPERTS: support@blirt.eu

		Agarose gel	Bacterial culture	Blood	Buffy coats	Cell culture	Enzymatic reaction mixture	Formalin/Paraffin embedded tissue	Hair / Insect	Physiological fluids	Plasma / Serum	Plasmid vectors	Solid tissue	Semen	Swabs – buccal, nasal, pharyngeal, vaginal	Yeast	Homogenisation bead beating tubes	Endotoxin-free	low elution volume / micro spin columns	DNase I present	Concentrator
GENOMIC DNA Isolation Kits																					
EM13	GENOMIC DNA KIT – universal		✓	✓	✓	✓		✓	✓	✓	✓		✓	✓	✓	✓					
EM02	DNA BACTERIA KIT		✓																		
EM03	DNA TISSUE KIT					✓		✓	✓	✓			✓								
EM04	DNA TISSUE PLUS KIT					✓		✓	✓	✓			✓				✓				
EM05	DNA BLOOD KIT			✓	✓						✓	✓									
EM06	DNA SWAB & SEMEN KIT													✓	✓						
RNA Isolation Kits																					
EM09.1	TOTAL RNA KIT					✓							✓								✓
EM11.1	TOTAL RNA PLUS KIT					✓							✓				✓				✓
EM12	miRNA KIT					✓							✓								
EM31	TOTAL RNA MICRO SPIN KIT					✓							✓						✓	✓	✓
EM15	RNA & DNA KIT					✓							✓								
EM25	RNA BACTERIA & YEAST KIT		✓												✓						✓
EM30	EXTRAZOL		✓			✓							✓								
PLASMID DNA Isolation Kits																					
EM01.1	PLASMID MINI KIT		✓										✓								
EM16/17	PLASMID MIDI KIT		✓										✓					✓			
EM18/19	PLASMID MAXI KIT		✓										✓					✓			
DNA FRAGMENTS Purification Kits																					
EM07.1	DNA CLEAN-UP KIT						✓														
EM08.1	DNA GEL-OUT KIT	✓					✓														
EM26.1	DNA CLEAN-UP & GEL-OUT KIT	✓					✓														
EM28	DNA CLEAN-UP & GEL-OUT MICRO SPIN KIT	✓					✓												✓		✓